

Ambu® Baby

40

The robust Ambu Baby is a lifelike simulation of a baby up to the age of one year. It is the youngest member of the Ambu training family and has all the sought-after Ambu qualities but with additional features to provide efficient training:

Ambu Baby Features

- The head can be both rotated and extended but the airway will only open in the correct (sniff) position
- Stomach ventilation can be visually checked

- Realistic anatomy allows the trainee to find the correct compression point, making the chest/abdomen rise realistically
- Realistic brachial pulse provides correct training
- Obstruction of the airway by an on/off slider can be simulated by the instructor
- The Ambu patented hygienic system facilitates knowledge retention through the look-listen-feel concept whilst providing trainee comfort and safety



ORDERING INFORMATION

CAT. NO.

Ambu Baby Manikin, for CPR with 5 face pieces and 100 head bags supplied in a softcase	256 001 000
Head bag (supplied in sets of 100)	256 000 702
Face piece (supplied in sets of 5)	256 000 703

Ambu® Man

41

The Ambu Man is an adult-sized training manikin that provides an exceptionally lifelike representation of the human anatomy, particularly those features important to training in modern resuscitation.

Ambu Man Features

- Ambu's unique hygienic system protects the students from cross-infection, makes internal cleaning entirely unnecessary and facilitates the look-listen-feel concept
- The stiffness of the chest is infinitely variable, allowing students to train on different body builds
- Airways open only when head is correctly hyper-extended
- A carotid pulse can be activated manually by the instructor

- A carotid pulse is also activated during correct chest compressions
- The mechanical monitoring instrument gives you instant feedback on ventilation volume, stomach inflation, chest compression depth in mm and wrong hand position
- PC connection is possible, allowing connection of Ambu CPR Software Kit, Ambu CPR Printer or the new USB Interface (see page 24)
- An upgrade kit can be purchased so that your BLS manikin can also function as an early defibrillation manikin
- Easy change of scale plates on the mechanical monitoring instrument if guidelines change



ORDERING INFORMATION

CAT. NO.

AmbuMan Model C, Torso with hooks for arms, computer interface, 100 head bags and 5 face pieces, supplied in jogging suit jacket and training mat.	234 003 000
AmbuMan Model CA, Torso with arms, computer interface, 100 head bags and 5 face pieces, supplied in jogging suit jacket and training mat.	234 023 000
AmbuMan Model CF, Fullbody, computer interface, 100 head bags and 5 face pieces, supplied in jogging suit jacket and training mat.	234 013 000

