

Ambu® Baby R resuscitator

34

The Baby R is a latex-free, re-usable resuscitator designed for the manual ventilation of babies, infants and children with a max. body weight of approx. 20 kg (4-5 years).

The suitably small dimensions of the Baby R means it can be used in incubators and can easily be operated with one hand. The double walling not only maximizes

user comfort by reducing hand fatigue but also imitates lung resistance.

Because the Baby R is fully re-usable it can be repeatedly autoclaved at 134° C (except the oxygen reservoir and the face mask). It can deliver a maximum tidal volume of 300 ml but it comes with a pressure limitation valve for patient safety. The volume of oxygen that can be held in the reservoir tubing is approx. 100 ml.



TECHNICAL DATA

Max. tidal volume of approx 300 ml
Patient connector: 22/15 mm (ISO)
Weight without mask: 160 g

Size: (Length x Diameter) 240 x 75 mm
Volume O₂ Reservoir: approx. 100 ml
Expiratory connector: 18 mm

ORDERING INFORMATION

	CAT. NO
Baby R with Paedi Valve and O ₂ Reservoir tube	083 025 000
Baby R with Paedi Valve, O ₂ Reservoir tube and face mask OA	083 019 000
Baby R with Paedi Valve, silicone transparent face blue cuff mask size 0 and O ₂ Reservoir tube	083 021 000
Baby R with Paedi valve with pressure limitation, O ₂ Reservoir tube and transparent face mask size OA	083 023 000



Ambu® Disposable Face Masks

Cost Effective And Efficient

The Ambu Disposable Face Masks are designed for use with manual and automatic ventilators and resuscitators.

The face masks have a thumb rest for user comfort and a soft, shaped, air-filled cuff to facilitate a tight fit to a patients face with an easy grip for the rescuer.

Because of the single patient use concept, there is no need for time consuming cleaning, sterilizing and disinfecting procedures. The added risks of cross infection is eliminated.

All versions, apart from the neonate and infant sizes, are supplied with a coloured hook ring so that the mask can be prepared for fastening by a head strap. The hook ring can easily be removed if it is not required.

COMPLETELY DISPOSABLE!

ORDERING INFORMATION

	CAT. NO.
Large Adult Disposable Face mask	000 252 056
Medium Adult Disposable Face mask	000 252 055
Small Adult/Child Disposable Face mask	000 252 054
Toddler Disposable Face mask	000 252 053
Infant Disposable Face mask	000 252 052
Neonate Disposable Face mask	000 252 051

