

# Ambu® Cardiac Care Trainer System

44

The Ambu Cardiac Care Training Manikin has been designed for training Advanced Life Support in ACLS algorithms and teamwork, either indoors or on location in field situations.

## Features

- Ambu Cardiac Care is modular expandable, meaning that you can invest gradually, updating the system as you go along
- The Manikin has the well tested mechanical monitoring system that gives you instant feedback on ventilation volume, stomach inflation, chest compression depth in mm and incorrect hand positioning
- The chest rigidity is adjustable so that different body builds can be simulated
- Ambu Cardiac Care has flexible defib electrodes so that you can train on this manikin in conjunction with with all leading brands of defibrillators

- Intubation with ordinary oral tube, laryngeal mask and Combitube is possible
- 4 ECG leads
- If the IV trainer (see page 45) is attached, diversion and IV training can be taught and practiced, making the whole scenario more realistic
- Ambu CPR Software Kit and Printer or the New USB interface can be connected

## Ambu ECG Box

- Ambu ECG Box simulates as many as 26 different rhythms and 2 artefacts, giving a wide range of ECG options
- The ECG Box has two optional power sources: an internal battery pack or an external 9V AC adapter giving you the flexibility to train on location and the security of knowing that the system will always work even if the batteries are low or missing



Ambu ECG Box has two optional power sources (internal battery or external 9V adapter)



## ORDERING INFORMATION

## CAT. NO.

Ambu Cardiac Care Trainer System, fullbody incl. jogging suit  
ECG box and hard case

284 013 000

Ambu ECG Box inc. Adapter and 7 x 1,5V batteries

284 000 730