

# Ambu® CPR Pal®



## Application

The low-priced Ambu CPR Pal is a simple but effective training manikin ideal for ACD CPR with the Ambu CardioPump®.

## Features

- The Ambu patented hygienic system protects the student and makes internal cleaning entirely unnecessary.
- Look – Listen – Feel. Just as in real life, ventilation can be seen and exhalation through the nose and mouth can be felt and heard.
- Lightweight and sturdy, hence very easy to transport.
- Airways open only when head is correctly hyperextended.
- Stomach inflation indicated by a whistling sound.
- ACD CPR-compatible (optional ACD CPR base board).



Ambu Hygienic System

## Specifications

- Manikin weighs 3.8 kg with and 2.7 kg without carry case.
- Total manikin length is 73 cm.

Ambu CPR Pal ordering specifications	Cat. No.
---	----------

Ambu CPR Pal with 5 mouth/nose pieces, 50 head bags and 1 T-shirt supplied in a soft carry case	259 001 000
Set of 100 head bags	000 071 017
Set of 25 mouth/nose pieces	071 000 586
ACD CPR base board for CPR Pal	259 000 701

For further details please contact your local Ambu representative.