

Ambu® LifeKey

The CPR Barrier You Can Take Everywhere

The Ambu LifeKey is a proven, easy to use, and inexpensive barrier with a one-way valve that prevents direct contact with a patient's mouth, nose and face and helps overcome any hesitations to start resuscitation.

The Ambu LifeKey helps protect the rescuer and increases bystander CPR in emergency situations. The face shield and one-way valve is semitransparent allowing the rescuer to visually check the patient's saliva or blood. Standard resuscitation techniques can be used without modification.

The small size of the Lifekey means that it can always be carried around, on a key ring, in a pocket or in a first aid kit. It is suitable for any patient, any age group. An elastic ear band fits over the patient's ears to keep the mask in place.

Folded it measures approx. 30 x 20 x 10 mm. It is available in a hard or soft case.

The Ambu Holster Pack is a soft case with belt loops that contains the disposable, easy to use face shield and disposable gloves.



ORDERING INFORMATION

	CAT. NO.
LifeKey in black softcase	248 206 102
LifeKey in blue softcase	248 206 103
LifeKey in red softcase	248 206 104
LifeKey in yellow softcase	248 206 105
LifeKey in green softcase	248 206 106
LifeKey in orange softcase	248 206 108
LifeKey in white hard case	248 204 101
LifeKey in red hard case	248 204 104
Holster Pack containing LifeKey, one set of gloves alcohol pad	248 510 102

Ambu® LifeKey private label

TECHNICAL DATA

Minimum Order Quantity: 250

Printing area (logo size) hard case: 20 x 20 mm

Available colours: black, blue, green, red, white and yellow

Printing area (logo size) soft case: 14 x 39 mm

Available colours: black, blue, red, yellow, green and orange and white hard case

ORDERING INFORMATION

	CAT. NO.
Ambu LifeKey in soft case Private label X=2 (black), X=3 (blue), X=4 (red), X=5 (yellow), X=6 (green), X=8 (orange)	248 206 30X
Ambu LifeKey in hard case Private label X= 1 (white), X=2 (black), X=3 (blue), X=4 (red), X=5 (yellow), X=6 (green), X=7 (grey)	248 204 30X

