

# Ambu® Military Mark III



## Designed For Toxic Atmospheres

The Military Mark III is similar in design to the standard Mark III (see page 9) but it has been further developed to comply with military specifications.

- It has a butyl outer cover for optimal resistance to toxins during ventilation in volatile atmospheres
- It has a gas filter adapter suited to all known gas filters to provide ventilation with filtered air and supplemental oxygen

## ORDERING INFORMATION

## CAT. NO.

|  |             |
|--|-------------|
| Military Mark III resuscitator complete with butyl outer cover, patient valve, gas filter adapter, (DIN norm 3182) adult transparent face mask, and special hanging strap with snap hook | 233 003 000 |
|--|-------------|

|  |             |
|--|-------------|
| Military Mark III resuscitator complete with butyl outer cover, patient valve, gas filter adapter, (US version, filter thread 2.3 in x 8 threads/in) adult transparent face mask, and special hanging strap with snap hook | 233 006 000 |
|--|-------------|

# Ambu® Transparent Face Masks



## 100% Natural Rubber Latex

The Ambu Transparent Face Masks, similar in design to the Transparent Silicone Face Masks, are designed for use with manual and automatic resuscitators and ventilators.

Since the Face masks are made from 100% natural rubber latex, they are the ideal solution for manual ventilation in toxic atmospheres, e.g fire and chemical accidents.

The masks are re-usable and can therefore be autoclaved repeatedly at 121° C.

## ORDERING INFORMATION

## CAT. NO.

|   |             |
|---|-------------|
| Adult Ambu Transparent Face Mask, size 5  | 000 017 000 |
| Adult/large child Ambu Transparent Face Mask, size 4                                  | 000 012 000 |
| Child Ambu Transparent Face Mask, size 2  | 000 016 000 |
| Small children Ambu Transparent Face Mask, size 0                                     | 000 015 000 |
| Small children Ambu Transparent Face Mask, size 0 with boring                         | 000 014 000 |
| Infant Ambu Transparent Face Mask, size 0A (cuff is made from chloroprene latex free) | 000 013 000 |