

Ambu® Open Cuff Silicone Face Masks

100% Latex Free

The latex-free Ambu Re-usable Open-Cuff Silicone Face Masks are designed for use with manual and automatic resuscitators and ventilators.

The face masks are transparent to allow for visual

check of bleeding, vomitus and spontaneous breathing.

These face masks are the ideal solution where economy is a priority.

The masks can be repeatedly autoclaved at 134° C.



ORDERING INFORMATION

CAT. NO

Adult Open Cuff Silicone Face Mask, size 5	000 251 006
Adult/Large Child Open Cuff Silicone Face Mask, size 4	000 251 005
Child Open Cuff Silicone Face Mask, size 3	000 251 004
Child Open Cuff Silicone Face Mask, size 2	000 251 003
Infant Open Cuff Silicone Face Mask, size 1	000 251 002
Infant Open Cuff Silicone Face Mask, size 0	000 251 001

Ambu® Transparent Silicone Face Masks

100% Latex Free

Long Lasting And Self-Inflating

The Ambu Transparent Silicone Face Masks have been designed for use with manual and automatic resuscitators and ventilators.

The Silicone Face Masks are self-inflating and the cuff pressure can be regulated for an easy and efficient seal. They are comfortable to grip and feature a thumb rest so rescuers can achieve a tight seal. Visual checks of bleeding, vomitus and

spontaneous breathing can be made through the transparent dome.

The Ambu Transparent Silicone Face Mask, size 0, with boring, impedes ventilation pressures higher than approx. 40 cm H₂O, enabling the rescuer to administrate correct ventilation on small children with a body weight of at least 15 kg (3 years).

Because the Ambu Transparent Silicone Face Masks are long lasting and re-usable they can be autoclaved repeatedly at 134°C.



ORDERING INFORMATION

CAT. NO.

Adult Silicone Transparent Face Mask, size 6	000 319 000
Adult Silicone Transparent Face Mask, size 5	000 317 000
Adult/large child Silicone Transparent face mask, size 4/3	000 312 000
Child Silicone Transparent face mask, size 2	000 316 000
Infant Silicone Transparent face mask, size 0	000 315 000
Infant Silicone Transparent face mask, size 0, with boring	000 314 000
New Transparent face mask, size 0a	000 013 000

See also page 18 for transparent face masks

