

Ambu® Disposable PEEP Valves



Precision And Stability

Ambu Disposable PEEP (Positive End Expiratory Pressure) Valves have been designed to provide unique resistance characteristics when used with manual resuscitators, ventilators, anaesthesia machines, (where specified by the manufacturer) and CPAP systems.

Ambu PEEP Valves can be fitted directly to the patient valve of any original Ambu resuscitator.

Because the PEEP Valves are disposable cross contamination issues are avoided.

Ambu PEEP Valve 20 is adjustable between 1.5-20 cm H₂O and can be supplied with an adapter with an O.D of 22 mm making them suitable for use with neonate resuscitators

TECHNICAL DATA

Size (Length x Diameter): 65 x 45 mm
Connection Size: I.D 30 mm

Weight approx.: 40 g
Adapter size: O.D 22 mm

ORDERING INFORMATION

CAT. NO.

Ambu Disposable PEEP Valve 20	199 102 001
Ambu Disposable PEEP Valve 20 with adapter	199 103 001

Ambu® Re-usable PEEP Valves

37



Precision And Stability

Ambu PEEP (Positive End Expiratory Pressure) Valves have been designed to provide unique resistance characteristics when used with manual resuscitators, ventilators, anaesthesia machines, (where specified by the manufacturer) and CPAP systems.

Ambu PEEP Valves can be fitted directly to the patient valve of any original Ambu resuscitator without additional connectors.

Ambu PEEP Valve 10 is adjustable between 0-10 cm H₂O and Ambu PEEP Valve 20 is adjustable between 1.5-20 cm H₂O.

The re-usable Ambu PEEP valves can be repeatedly autoclaved at 134°C.

TECHNICAL DATA

Size (Length x Diameter): 65 x 45 mm
Connection: 18/30 mm for resuscitators, 19 mm for anesthesia machines

Weight approx.: 40 g

ORDERING INFORMATION

CAT. NO

PEEP Valve 10, adults 30 mm and baby 18 mm connector	137 001 000
PEEP Valve 10, adult 30 mm	000 137 000
PEEP Valve 10, baby 18 mm	000 138 000
PEEP Valve 20, with inlet connector 30 mm	000 213 000
PEEP Valve 20, for respirators/CPAP system 22/19 mm	000 177 000