

Ambu® Perfit™ Extrication Collars

The Perfect Fit For Neutral Alignment

Ambu offers two variations of collar.

The **Ambu Perfit** and the **Ambu Perfit ACE** (Adjustable Collar for Extrication). Both types of collar are designed to assist with the maintenance of neutral alignment, prevention of lateral sway and anterior-posterior flexion and extension of the cervical spine during transport and routine patient care or movement.

All Ambu Extrication Collars have the following features

- Ambu collars come with a comprehensive sizing system for easy and accurate sizing
- One-piece extrication collars simplify application and at the same time there are no extra pieces that can get lost
- Nasal cannula holders are located on the collar body to secure the oxygen tubes (except for infant)
- All collars are radiotranslucent, CAT & MRI compatible and come with a ventilated posterior shell for fluid drainage and a large hole for easy access. All of these features allow for patient examinations to take place safely without removing the collar
- Safety buttons allow for maximum reliability
- Ambu collars are fully disposable, you no longer have to worry about collars being returned to the source of origin and cross contamination issues are avoided
- Ambu collars come with a pre-shaped flip-chin piece to ensure they can be stored completely flat. Furthermore the flip-chin piece facilitates intubation and removal of vomit from a patient's mouth



Ambu® Perfit

Ambu Perfit

The Ambu Perfit extrication collar is a one-size collar,

which comes in four adult sizes, a paediatric size and an infant size.



TECHNICAL DATA

Shell material: polyethylene

Sizes:

Colour: White/light grey

1 (infant), 2 (paediatric), 3 (neckless), 4 (short), 5 (regular) and 6 (tall)

ORDERING INFORMATION

CAT. NO.

Perfit Infant: size 1	264 201 000
Perfit Paediatric: size 2	264 202 000
Perfit Neckless : size 3	264 203 000
Perfit Short: size 4	264 204 000
Perfit Regular: size 5	264 205 000
Perfit Tall: size 6	264 206 000
Carry bag for collars, fits 6 pieces	264 008 000
Sample pack, one of each Perfit size	264 020 000