

Ambu® Res-Cue Pump



Power In The Palm Of Your Hand

A cost effective alternative to electrical and foot powered devices, the Ambu Res-Cue Pump is a portable, hand-operated emergency suction pump.

The Ambu Res-Cue Pump has an adjustable vacuum which makes it possibly the most versatile manual suction pump available and it is probably the safest and easiest to use too.

The Overfill-Protection System prevents contamination of the main unit and averts any inconvenient

splashing of fluids. It is very easy to squeeze, requiring little force to obtain effective suction.

This lightweight yet robust pump utilises a unique disposable collection container which eliminates the need for cleaning and maintenance. Each disposable collection container is supplied with two catheters for adults and paediatrics, ensuring you are prepared for any eventuality.



TECHNICAL DATA

Vacuum: 450 mm Hg
 Max. flow: > 20 l/min.
 Container capacity: 300 ml (disposable)
 Supplied with:

Adjustable vacuum: (100 % or 50%)
 Size: 206 x 96 x 104 mm
 Weight: 260 g

Two suction tubes, 6 and 13 mm and 2 catheter-adapters, one for finger-tip connection and one for funnel shape connection

NEW!
 Now with two
 catheter-adapters



ORDERING INFORMATION

CAT. NO.

ResCue Pump including, two suction tubes (6 and 13 mm) and two catheter adapters	276 000 001
Disposable collection container	276 000 010
ResCue Pump Carry Case	276 000 011

Ambu® Suction Booster

Additional reservoir

The Ambu Suction Booster significantly increases the speed of suction and makes it possible to remove larger solid particles which could obstruct average suction tips.

It is easily connected to the suction tip of all Ambu Suction Pumps.

Ambu Suction Booster has a reservoir capacity of 250ml and is autoclavable at 121° C.



ORDERING INFORMATION

CAT. NO

Ambu Suction Booster	136 001 000
----------------------	-------------