

Ambu® ResCue Mask

Professional Protection

The Ambu ResCue Mask is designed to provide protection for the rescuer when manually resuscitating patients. It can also be used with a manual resuscitator.

Ambu ResCue Mask has a one-way valve with bacterial filter that helps to protect the rescuers in an emergency situation. It prevents direct contact with a patient's mouth, nose and face and helps overcome any hesitations to start resuscitation.

The transparent dome allows the rescuer to visually check the patient's airway status and the oxygen inlet (optional) facilitates increased oxygen supply to the patient. To help keep the mask in place the mask is provided with an elastic band to place around the patient's head.

It measures 13.8 x 11 x 4.8 cm. Available in hard case, soft case or plastic bag.

ORDERING INFORMATION

CAT. NO.

Ambu Res-Cue Mask - Hard red case with O ₂ inlet	000 252 225
Ambu Res-Cue Mask - Hard red case without O ₂ inlet	000 252 226
Ambu Res-Cue Mask - Soft case with O ₂ inlet	000 252 223
Ambu Res-Cue Mask - Soft case without O ₂ inlet	000 252 224
Ambu Res-Cue Mask - In plastic bag with O ₂ inlet	000 252 221
Ambu Res-Cue Mask - In plastic bag without O ₂ inlet	000 252 222
One way valve with bacterial filter	000 252 253
Ambu Res-Cue Mask - Adult and Infant mask in soft case without O ₂ inlet	000 252 123



Ambu® ResCue Mask private label

TECHNICAL DATA

Organisation, company or sponsor logo

Logo Size: Max. 50 x 40 mm

Colour: Four colour print is available on either a white or transparent sticker
1 colour print, if logo is printed on the case itself

The hard case is available in all pantone colours

Min Order: 250 pieces

ORDERING INFORMATION

CAT. NO.

Ambu ResCue Mask Private label - Hard Case, without O ₂ inlet	000 252 236
Ambu ResCue Mask Private label - Hard Case, with O ₂ inlet	000 252 235

