

# Accessories & Spare Parts

## Ambu® Mark III Resuscitator

ACCESSORIES	CAT. NO.
1 Tube adapter, 24 mm female for extension tube	209 000 507
2 Tube adapter 28/32 mm, 28 mm female connector for fitting on inlet connector of Mark III patient valve, 32 mm male connector for fitting the inlet valve of Mark III bag	209 000 506
3 Transparent face mask	see page 9
4 Ambu Single Shutter Patient Valve w/swivel, 24 mm inlet connector, 22/15 mm patient connector and 30 mm outlet connector. Fits both Ambu Mark III and Ambu Silicone resuscitators	295 000 501
5 Gas filter adapter	211 000 501
6 Mark III patient valve	000 208 000
7 Oxygen reservoir bag	210 000 701
8 PEEP valve 20, for adult	000 213 000
9 PEEP valve 10, for adult resuscitators with connector, ID 30 mm, for adult patient valves	000 137 000
10 Spiral hose for O <sub>2</sub> Reservoir, EPDM rubber	209 000 508
11 Carry plastic bag	209 000 707
12 Ambu Pack with transparent cover	230 000 001
14 Extension tube for Mark III, 1068 mm	209 000 703
15 Extension tube for Mark III, 305 mm	209 000 701
16 Extension tube for Mark III, 915 mm	209 000 702

SPARE PARTS	CAT. NO.
17 Hanging strap with identification tag	209 000 501
18 Special hanging strap with snap lock	215 000 501
19 Outer cover	209 000 502
21 Butyl outer cover	233 000 502
22 Inlet valve housing, complete	209 000 510
24 Inlet valve membrane	209 000 501
27 Inner bag, complete	209 000 504

