

Ambu® Silicone Resuscitator

ACCESSORIES

	CAT. NO.
1 Tube adapter, 24 mm female for extension tube	209 000 507
2 Tube adapter 28/32 mm, 28 mm female connector for fitting on inlet connector of Mark III patient valve, 32 mm male connector for fitting the inlet valve of Mark III bag.	209 000 506
3 Oxygen reservoir bag	210 000 701
6 Ambu Single Shutter Patient Valve w/swivel, 24 mm inlet connector, 22/15 mm patient connector and 30 mm outlet connector. Fits both Ambu Mark III and Ambu Silicone resuscitators	295 000 501
7 Oxygen reservoir bag, reusable	245 000 519
8 Oxygen reservoir tube, 250 mm	288 000 510
9 Patient valve complete, neonate	288 000 501
10 PEEP valve 20, for adult	000 213 000
11 PEEP valve 10, for adult resuscitators with connector, ID 30 mm, for adult patient valves	000 137 000
12 PEEP valve 10, for baby resuscitator, 18 mm	000 138 000
13 Extension tube for Mark III, 1068 mm	209 000 703
14 Extension tube for Mark III, 305 mm	209 000 701
15 Extension tube for Mark III, 915 mm	209 000 702

SPARE PARTS

	CAT. NO.
16 Inlet valve housing with valve disc and hanging strap (grey)	245 000 506
18 Flange nut (grey)	245 000 507
20 Outer flange for inlet valve	245 000 516
22 Silicone resuscitator bag, infant/child complete without patient valve and mask	246 000 504
23 Silicone resuscitator bag adult with inlet valve	245 000 502
24 Resuscitator bag neonate, complete without patient valve and mask	288 000 504
25 Resuscitator bag neonate with support ring	288 000 503
26 Resuscitator bag infant/child with inlet valve	246 000 502
28 Inlet valve cover	288 000 507
29 Inlet valve housing with three valve shutters (adult)	245 000 501
35 Valve membrane (neonate) 25 min.	288 000 506
37 Valve disc for inlet valve (adult) 25 min.	245 000 514

