

# Ambu® TwinPump

SPARE PARTS	CAT. NO.
1 Suction tube, 1350 mm with combination suction tip and angle Connector	239 000 509
2 Pump chassis, complete with valve diaphragms incl. O-rings and piston rings. 1000 ml container	240 000 501
3 Pump chassis, complete with valve diaphragms incl. O-rings and piston rings. 600 ml container	239 000 501
4 Bottom cover	239 000 514
5 Pedal with retaining springs	239 000 508
6 Piston ring	239 000 516
7 O-ring	239 000 518
8 Combination suction tip	239 000 510
9 Valve diaphragm, separate 10 pcs. min.	239 000 519
10 Manifold pipe	239 000 517
11 Angle connector	239 000 515
12 Cylinder with drawlink and valve diaphragm	239 000 504

